Al-based Inline Prediction of Tablet Properties for Accelerating Formulation Development

Nico Kühnel¹, Jan Streckenbach¹, Pierre Braun¹, Moritz Rosch¹, Friederike Gütter¹, David Dormagen^{2,3}, Tim Landgraf^{2,3}.

¹KORSCH AG, Breitenbachstraße 1, 13509 Berlin, Germany

² CIATA UG (haftungsbeschränkt), Arnimallee 7, 14195 Berlin, Germany

³ Freie Universität Berlin, Arnimallee 7, 14195 Berlin, Germany

INTRODUCTION AND AIM

critical quality attributes (CQAs) such as tablet hardness, tensile strength, process features during compaction. It enables faster screenings of blends solid fraction, and tablet weight. These CQAs are traditionally obtained by and process settings while reducing reliance on destructive end-testing. destructive offline tests that consume material and slow down iteration.

Development of solid oral dosage forms depends on timely knowledge of A machine learning model was trained to predict these properties inline from

MATERIALS AND METHODS

The model development followed a structured data science workflow¹. First, targets (CQAs). This was followed by data understanding, including exploratory analysis of datasets with more than 300,000 tableting cycles, 200 punch profiles and 30 machine simulation profiles. During this step, erroneous, implausible or inconsistent records e.g., missing reference values or corrupted sensor streams were identified and excluded. In the data preparation phase, relevant force-displacement and tablet energy features like plastic, elastic and ejection energy were extracted, engineered and normalized. Subsequently, different machine learning regression models were trained and benchmarked in the modelling phase. A model performance of R² with 0.95 and mean absolute percentage error of 15% was achieved using

time-separated splits to avoid data leakage. The model was validated with business and process understanding were established to define prediction an external dataset of 9,364 tablets from 73 blends, using 36 machine and 67 punch profiles.

> In a representative case study, an API-containing customer formulation was compressed on a STYL'One Evo equipped with 12 mm round EU-B tooling and configured to emulate a KORSCH XL 400⁴ at 100 min⁻¹ turret speed. 12 tablets were compressed at a main compaction force ranging from 6.7 to 13.5 kN. The machine learning model was used to predict the CQAs directly from compaction data. Destructive reference measurements (hardness, weight, diameter, height) were taken on a Kraemer LAB.line P5 tablet tester and compared to the predictions.

> > --- ideal (measured = prediced)

hardness average ± standard deviation

± 15 % range

RESULTS & DISCUSSION

Using a prediction model reduces reliance on destructive tests and provides actionable, early-stage feedback on hardness, tensile strength, solid fraction, and weight. The relative deviation of predicted hardness vs. measured hardness across the validation set is centered around zero percent, but reveals increased errors at the boundaries probably caused by training data sparsity in extreme regions and historical unsystematic acquisition patterns (Figure 1).

The predicted vs. measured scatterplots for an API-containing customer product confirm high agreement in the investigated range (Figures 2 and 3), while the graphical user interface illustrates immediate usability in development workflows (Figure 4).





Figure 1: Relative deviation of predicted hardness vs. measured hardness for validation set of 9,364 tablets.



§ 700 \vdash 660 640 620

Figure 3: Predicted vs. measured weight for an API-containing customer product.

measured weight in mg







Figure 2: Predicted vs. measured hardness for an API-containing customer product.

--- ideal (measurment = predicted)

weight average ± standard deviation

± 10 % range

760



OUTLOOK

Future work will expand and refine the machine learning model prediction by densifying training data in boundary regions, where current predictions are less accurate. Furthermore, adopting an ensemble composition of neural networks is expected to improve robustness². Finally, the CQAs will be expanded to include disintegration time and friability, providing a more comprehensive prediction framework for pharmaceutical development.